

Ancillary component designed to cushion shock and shear forces caused by vertical and rotational movement. Ideal for sports and activities such as golf, tennis and hiking.

TTPro

## Specifications

<b>Max. User Weight/Impact Level:</b>	100kg/ 1-4 125kg/ 1-3
<b>Component Weight:</b>	From 317g
<b>Build Height:</b>	135mm - 260mm (pyr) 145mm (clamp)
<b>System diameter:</b>	30mm



Pyramid



Clamp



Suitable for outdoor use

## Order Example


Pyramid			
Product Code	Rod	Spring Set	
<b>TTPRO</b>	<b>PYR</b>	<b>3</b>	<b>B</b>

eg. TTPRO-PYR3B

Clamp			
Product Code	Rod	Spring Set	
<b>TTPRO</b>	<b>30</b>	<b>3</b>	<b>B</b>

eg. TTPRO-30-3B

## Build Height



Impact Level	User Weight								Rod and Spring Set
	44-52 100-115	53-59 116-130	60-68 131-150	69-77 151-170	78-88 171-195	89-100 196-220	101-116 221-255	117-125 256-275	
Low	1P	2P	2P	2W	2W	3B	3B	3B	
Moderate	2P	2W	3W	3W	3B	3B	3B	3B	
High	3W	3W	3B	3B	3B	3B	-	-	

41  
32