

## DESCRIPTION

A dynamic response foot for activity level 2-3 users and designed for use with the Multiflex Ankle, the DR2 has a user adjustable carbon fibre toe spring that can be set for different terrain or activities.

## FEATURES

- User adjustable response
- Multi-axial function
- 20mm ±15mm heel height range fixed at limb assembly
- Full length toe lever




## SPECIFICATIONS

Max. user weight:	125kg
Activity level:	2-3
Size range:	22-30
Component weight:	430g†
Build height:	55mm
Heel height:	20 +/-15mm

## PART NUMBERS

Size	Left	Right
22cm	529120	529121
23cm	529122	529123
24cm	529124	529125
25cm	529126	529127
26cm	529128	529129
27cm	529130	529131
28cm	529132	529133
29cm	529134	529135
30cm	529136	529137

## BUILD HEIGHT



### Multiflex Ankle Standard 30mm and 35mm

	30mm	35mm
Ankle Kits:	379342 (W)	409125 (W)
	379343 (O)	409126 (O)
	379344 (B)	409127 (B)
Snubber Kits:	409007 (W)	
	409008 (O)	
	409009 (B)	

### Multiflex Ankle Pyramid

Ankle Kits:	409353 (W)
	409354 (O)
	409355 (B)
Snubber Kits:	409007 (W)
	409008 (O)
	409009 (B)

DR2 -  
Multiflex ankle and snubber selection.

Activity	User Weight							kg lbs	Ball and snubber
	44-52 100-115	53-59 116-130	60-68 131-150	69-77 151-170	78-88 171-195	89-100 196-220	101-125 221-275		
2	WW	WW	WW	WO	OO	OB	BB		
3	WW	WW	WO	OO	OB	BB	BB		

†Component weight shown is for a size 26cm