

Endolite Comfort Liner - Cushion



Features

- Trans-tibial
- Seamless bi-directional stretch to improve comfort and contouring
- Tendresse™ finish reduces shear stress on the skin
- 6 month warranty

Part Numbers

CTTCPXXL Comfort - Cushion

XX indicates size, see sizing chart and illustrations for details

Sizing Chart

Measurements taken 4cm from the distal end

Measured Size (cm)	Part Number	Description
22cm	CTTCP22L	COMFORT - CUSHION 22
23.5cm	CTTCP23L	COMFORT - CUSHION 23.5
25cm	CTTCP25L	COMFORT - CUSHION 25
26.5cm	CTTCP26L	COMFORT - CUSHION 26.5
28cm	CTTCP28L	COMFORT - CUSHION 28
30cm	CTTCP30L	COMFORT - CUSHION 30
32cm	CTTCP32L	COMFORT - CUSHION 32
34cm	CTTCP34L	COMFORT - CUSHION 34
36cm	CTTCP36L	COMFORT - CUSHION 36
40cm	CTTCP40L	COMFORT - CUSHION 40

*note- 23.5 = CTTCP23L , 26.5 = CTTCP26L

Endolite Comfort Liner - Locking



Features

- Trans-tibial
- Pin lock version suitable for a variety of lock and pin preferences
- Stabilising external matrix with integrated stretch zones
- Tendresse™ finish reduces shear stress on the skin
- 6 month warranty

Part Numbers

CTTLPXXL Comfort - Locking

XX indicates size, see sizing chart and illustrations for details

Sizing Chart

Measurements taken 4cm from the distal end

Measured Size (cm)	Part Number	Description
22cm	CTTLP22L	COMFORT - TT LOCKING 22
23.5cm	CTTLP23L	COMFORT - TT LOCKING 23.5
25cm	CTTLP25L	COMFORT - TT LOCKING 25
26.5cm	CTTLP26L	COMFORT - TT LOCKING 26.5
28cm	CTTLP28L	COMFORT - TT LOCKING 28
30cm	CTTLP30L	COMFORT - TT LOCKING 30
32cm	CTTLP32L	COMFORT - TT LOCKING 32
34cm	CTTLP34L	COMFORT - TT LOCKING 34
36cm	CTTLP36L	COMFORT - TT LOCKING 36
40cm	CTTLP40L	COMFORT - TT LOCKING 40

*note- 23.5 = CTTLP23L , 26.5 = CTTLP26L